

Air Techniques, Inc. Introduces AirStar® and STS™ Enhancements

Date: 06/22/2009



MELVILLE, NEW YORK. Air Techniques, Inc. introduces the latest enhancements to its AirStar Air Compressors and STS Dry Vacuums. These new advancements were created to provide you with the best possible results for your practice.

In addition to the new dryer membrane technology that will soon be available on all AirStar models, we have worked to make the AirStar quieter than ever. To achieve minimal noise level AirStar Compressors now have improved intake filters, a redesigned purge muffler, new rubber shock mounts, and a flat frame to reduce vibration

The STS Dry Vacuum enhancements increasing product life, cleanliness, and space. A new horizontal motor design results in a smaller footprint for space-saving flexibility. Nickel-plating on the propeller and housing reduces oxidation therefore increasing efficiency and motor life. The increased size of the motor bearings ensures fewer failures and increased reliability. The CAS Washout Timer Control keeps debris from building up in the tank and helps to keep the tank clean by automatically sending water to the spray nozzle even when the pump is turned off.

For more information, visit the [Products page](#) of our website or call 1-800-AIR-TECH. Founded in 1962, Air Techniques is a leading manufacturer of quality equipment for the dental industry.
